


Laparoscopy

	300/5.0 + 300/10.0 Rigid Laparoscopes				
<p>Completely made of high quality steel for highest corrosion resistance. All mechanical bonds are laser welded. Distal and proximal cover lenses consist of sapphire glass and are soldered in. That means: no adhesive surfaces, no leakage and no scratching. "High-resolution" optic systems for special brilliant image quality. Althea endoscopes are compatible with shafts of all different manufacturers and can be integrated in your present equipment / accessories without any problems. Different lines of vision are marked in terms of color. The blue ring shows that the endoscope is autoclavable.</p>					
<p>>> Cover lens made of sapphire glass >> Corrosion-resistant >> 134°C autoclavable >> Brilliant image >> Absolutely leak proof</p>					
Diameter	5,00 mm	5,00 mm	10,00 mm	10,00 mm	10,00 mm
Working length	300 mm	300 mm	300 mm	300 mm	300 mm
Angle of view	0°	30°	0°	30°	45°
Aperture angle	75°	75°	75°	75°	75°
Article no.:	20 05 00	20 05 30	20 10 00	20 10 30	20 10 45