

Orthopedics



Completely made of high quality steel for highest corrosion resistance. All mechanical bonds are laser welded. Distal and proximal cover lenses consist of sapphire glass and are soldered in. That means: no adhesive surfaces, no leakage and no scratching. "High-resolution" optic systems for special brilliant image quality. Althea endoscopes are compatible with shafts of all different manufacturers and can be integrated in your present equipment / accessories without any problems. Different lines of vision are marked in terms of color. The blue ring shows that the endoscope is autoclavable.

- **Cover lens made of sapphire glass**
- **Corrosion-resistant**
- **134°C autoclavable**
- **Brilliant image**
- **Absolutely leak proof**

Compatible	STORZ	STORZ	STORZ	STORZ	Wolf	Wolf	Wolf
Diameter	4,00 mm	4,00 mm	4,00 mm	4,00 mm	4,00 mm	4,00 mm	4,00 mm
Working length	175 mm	175 mm	175 mm	175 mm	170 mm	170 mm	170 mm
Angle of view	0°	30°	45°	70°	0°	30°	70°
Aperture angle	95°	95°	95°	95°	95°	95°	95°
Article no:	10 40 00	10 40 30	10 40 45	10 40 70	13 40 00	13 40 30	13 40 70

Compatible	STORZ	STORZ	STORZ	STORZ	STORZ	STORZ	Wolf
Diameter	2,70 mm	2,70 mm	2,70 mm	2,70 mm	2,70 mm	2,70 mm	2,70 mm
Working length	187 mm	187 mm	175 mm	175 mm	110 mm	110 mm	180 mm
Angle of view	0°	30°	0°	30°	0°	30°	0°
Aperture angle	80°	80°	80°	80°	80°	80°	80°
Article no:	10 28 00	10 28 30	11 28 00	11 28 30	12 28 00	12 28 30	13 28 00

Althea Deutschland GmbH

Im Goldäcker 14 | D-88630 Pfullendorf | Tel. +49 (0) 7552 92801-0 | Fax. +49 (0) 7552 92801-11
E-Mail: de.info@althea-group.com | www.althea-group.com/de